



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Mettler-Toledo® InLab 731, conductivity electrode, with 1.2 m cable, mini DIN connector

€520.00
plus VAT &
Shipping

Product Images



Description

Conductivity sensor INLAB 731 - sensor for medium to high conductivities - very robust epoxy shaft - 4-pole graphite cell - cell constant: 0.57 cm⁻¹ - measuring range: 0.01...1000 µS/cm - connection: fixed cable with Mini-DIN plug - ideal for SevenMulti instruments - shaft length 120 mm - shaft diameter 12 mm

Additional Information

No.	TD-0457
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105004020
Transport temperature	Room temperature

