



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Mettler-Toledo® InLab 741, conductivity
electrode, with 1.2m cable, mini DIN connector**

**€735.00
plus VAT &
Shipping**

Product Images



Description

Conductivity sensor INLAB 741 – Sensor for low conductivities like clear and ultrapure water – very robust stainless steel shaft – 2-pole stainless steel (V4A) cell – cell constant: 0.105 cm⁻¹ – measuring range: 0.001....500 µS/cm – connection: Mini-DIN plug – ideal for SevenMulti instruments – shaft length 120 mm – shaft diameter 12 mm

Additional Information

No.	TD-0459
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105004024
Transport temperature	Room temperature

