



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Mettler-Toledo® InLab Cool, pH electrode

**€348.00**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

This combined pH electrode, whose reference electrode is filled with the electrolyte Friscolyte-B®, was developed for low-temperature samples (down to -30°C) and low-conductivity media. Specifications: Glass shaft (12 mm diam./ 120 mm long) S7 screw head Ceramic diaphragm ARGENTHAL™- reference with silver ion barrier pH range: 0...11 Temp. range: -30...80°C Incl. 25 ml reference electrolyte (FRISCOLYT-B®)

## Additional Information

---

No.	TD-0435
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105003174
Transport temperature	Room temperature

