



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Mettler-Toledo® InLab Expert Pt1000, pH electrode

€286.00
plus VAT &
Shipping

Product Images



Description

Ideal for a wide range of applications, including emulsions, suspensions, wastewater and solutions containing proteins and sulfides. Specifications: High quality plastic shaft (PEEK) and MultiPin™ head. Polymer electrolyte XEROLYT® Open connection ARGENTHAL™ - Reference Integrated temperature sensor (Pt1000) pH range: 0...14 Temp. range: 0...100°C

Additional Information

No.	TD-0424
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105003105
Transport temperature	Room temperature

