



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Mettler-Toledo® InLab Mono, pH half cell

€161.00
plus VAT &
Shipping

Product Images



Description

The use of a separate pH electrode (half cell) can save costs in those applications where the life of the reference electrode is much longer than that of the pH electrode. Specifications: Glass shaft (12 mm diam./ 120 mm long) S7 screw head pH range: 0...14 Temp. range: 0...100°C

Additional Information

No.	TD-0445
Manufacturer (Brand)	Mettler-Toledo
Transport temperature	Room temperature

