



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Mettler-Toledo® InLab Reach Pro-425, pH electrode

€411.00
plus VAT &
Shipping

Product Images



Description

Designed for accurate measurements on a wide range of different sample types. Now in an extra long version! Specifications: Glass shaft (12 mm diam./ 400 mm long) S7 screw head Ceramic diaphragm ARGENTHAL™-reference with silver ion barrier Integrated temperature sensor (NTC30 kΩ) pH range: 0...14 Temp. range: 0...100°C Incl. 25 ml reference electrolyte (3 mol/L KCl)

Additional Information

No.	TD-0416
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105003061
Transport temperature	Room temperature

