



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Mettler-Toledo® InLab Redox Ag

€351.00
plus VAT &
Shipping

Product Images



Description

Redox electrode consisting of a silver ring indicator electrode and a silver/silver chloride reference electrode. The measurement signal is generated on the surface of the precious metal by an electron exchange with the oxidation/reduction system of the measured medium. Specifications: Glass shaft (12 mm diam./ 120 mm long) S7 screw head ARGENTHAL™- reference with silver ion barrier Diaphragm: ceramic Metal: platinum ring Temp. range: 0...100°C

Additional Information

No.	TD-0452
Manufacturer (Brand)	Mettler-Toledo
Transport temperature	Room temperature

