



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Mettler-Toledo® InLab 738 ISM-10m, conductive electrode, with 10 m fixed cable, LTW connector

€720.00
plus VAT &
Shipping

Product Images



Description

InLab® 738 ISM-5m - as InLab 738 ISM # 51344110, but with 10 m fixed cable

Additional Information

No.	TD-0470
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105004114
Transport temperature	Room temperature

