



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Mettler-Toledo® InLab 742 ISM, conductivity electrode, with 2 m fixed cable, LTW connector

**€795.00**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

Conductivity sensor INLAB 742 - Sensor for low conductivities such as clear and ultrapure water - Very robust stainless steel shaft - 2-pole stainless steel (V4A) cell - Cell constant: 0.105 cm<sup>-1</sup> - Measuring range: 0.001....500 µS/cm - Connection: LTW connector - Ideal for SevenMulti instruments - Shaft length 120 mm - Shaft diameter 12 mm - Intelligent Sensor Management (ISM): sensor data, calibration, and maximum temperature are stored in the sensor and are automatically transferred to the meter when connected - Ideal for SevenGo and SevenCompact instruments

## Additional Information

---

No.	TD-0471
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105004116
Transport temperature	Room temperature

