



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Mettler-Toledo® InLab 742, conductivity electrode, with 5 m fixed cable, LTW connector

€790.00
plus VAT &
Shipping

Product Images



Description

Conductivity sensor INLAB 742 - Sensor for low conductivities like clear and ultrapure water - very robust stainless steel shaft - 2-pole stainless steel (V4A) cell - cell constant: 0.105 cm⁻¹ - measuring range: 0.001....500 µS/cm - connection: LTW plug - ideal for SevenMulti instruments - shaft length 120 mm - shaft diameter 12 mm

Additional Information

No.	TD-0477
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105004128
Transport temperature	Room temperature

