



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Mettler-Toledo® InLab 752-6mm

**€645.00**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

2-pole platinum micro conductivity sensor with 6 mm shaft diameter, integrated temperature sensor and fixed cable. Intended for use with METTLER TOLEDO's SevenCompact™, SevenMulti™ and FiveEasy™ bench-top meters, as well as the FiveGo™ portable meter for conductivity measurement of small sample quantities. Cell constant: 1.0 /cm. 2-pole design provides a uniform signal for accurate measurement of low to moderate conductivity samples. Glass body resists aggressive or organic sample types. Technical data: Shaft material glass Shaft length: 180 mm Shaft diameter 6 mm Temperature range °C: 0 ... 100 Measuring range 0.01 ... 112 mS/cm, cable and connector 1m, mini DIN plug

## Additional Information

---

No.	TD-0462
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105004031
Transport temperature	Room temperature

