



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Mettler-Toledo® InLab Expert Pro-ISM-10m, with 10 m fixed cable

€555.00
plus VAT &
Shipping

Product Images



Description

pH electrode InLab® Expert Pro-ISM for SevenGo instruments: - robust 3in1 electrode with PEEK shaft - polymer electrolyte - maintenance-free - contamination-resistant due to open connection - integrated temperature sensor - Intelligent Sensor Management (ISM): sensor data, calibration, and maximum temperature are stored in the sensor and are automatically transmitted to the meter when connected - fixed cable: 10m - connector: BNC (IP67), Cinch/RCA

Additional Information

No.	TD-0467
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105004104
Transport temperature	Room temperature

