



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Mettler-Toledo® InLab Flex-Micro, pH electrode

**€484.00**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

Combined micro pH electrode with 6mm shaft diameter and unique bendable plastic shaft => prevents breakage - ideal for small sample volumes and vessels, ... - Fixed cable with BNC connector - PTFE diaphragm - Argenthal™-reference system with silver ion barrier, Technical data: - pH range: 0 ... 14 - Temperature range °C: 0 ... 80 - Diaphragm resistance (25 °C): < 600 MOhm - Reference electrolyte: Friscoylt-B Bearing solution: 3M KCl - Shaft length: 180 mm - Shaft diameter: 6mm

## Additional Information

---

No.	TD-0433
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105003164
Transport temperature	Room temperature

