



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Mettler-Toledo® InLab Power Pro-ISM

**€660.00**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

InLab® Power Pro-ISM - combined pH electrode with unique SteadyForce™ reference system (DPA electrolyte with 2.5 bar overpressure = no contamination of the diaphragm possible) - contamination- and maintenance-free - ideal for viscous and liquid media such as dairy products, contaminated solutions, ... - integrated temperature sensor - MultiPin™ plug-in head - SteadyForce™- reference system with silver ion barrier - Intelligent Sensor Management (ISM): sensor data, calibration, and maximum temperature are stored in the sensor and are automatically transferred to the meter when connected - ideal for SevenGo and SevenCompact instruments, necessary connection cable: # 51344291 Technical data - pH range 0 ... 12 - Temperature range °C 0 ... 130 - Temperature sensor NTC 30 kOhm - Diaphragm resistance (25 °C) 600 MOhm - Reference electrolyte: DPA gel - Bearing solution: Fiscolyt B, #51350076 - Shaft length 170 mm - Shaft diameter 12 mm

## Additional Information

---

No.	TD-0481
Manufacturer (Brand)	Mettler-Toledo
VGKL number	105004211
Transport temperature	Room temperature

