



SI Analytics® Ids electrode Scppht-Sp-A-3M-Ids

Product Images

Description

IDS electrode glass shaft, insertion electrode, 3 m fixed cable with digital connector, 3 x ceramic diaphragm, electrolyte KCl 3 mol/l, Silamid® lead with silver ion barrier and double electrolyte, A-glass, length 120 mm, 12/8 mm Ø, -5...+100 °C, 0...14 pH

Additional Information

No.	SY-2435
Manufacturer (Brand)	SI Analytics

