



SI Analytics® Scienceline Plus Scppht-Mic-Amf-3M-Bnc-N

Product Images

Description

ScienceLine plus pH combination electrode Micro electrode, Pt 1000 temperature sensor, 3 m cable with BNC and banana plug, glass shaft, platinum diaphragm, electrolyte KCl 3 mol/l, cylindrical membrane, A-glass, length 70/130 mm, 12/5 mm Ø, -5...+10

Additional Information

No.	SY-2431
Manufacturer (Brand)	SI Analytics

