



Thermo Scientific Electrodes Orion Basic

Product Images

Description

Orion Basic electrodes Thermo Scientific Orion Basic electrodes are designed for Russell instruments and for simple pH measurements. Description: Refillable combination pH electrode with glass body For use with (instruments): Russell Meter Storage conditions: Store in Orion pH electrode storage solution 910001 Operating temperature (Anglo-American): 32 ° to 176 ° F Connection types: BNC connector Material: Glass shaft Operating temperature (metric): 0 °C° to 80 °C° Type: pH combination electrode, glass body pH range: 0 to 14 Unit Size: Each

Additional Information

No.	LD-2896
Manufacturer (Brand)	Thermo Scientific

