



Thermo Scientific Eutech Conductivity/ATC Electrodes

Product Images

Description

Eutech Conductivity/ATC Electrodes For fast and accurate conductivity measurements in the field, in laboratories and in schools. For use with (instruments): CON 700 / PC 700 / CON 2700 / PC 2700 Temperature range (metric): 0 °C to 80 deg;C Cell constant: K = 0.1 Type: electrode Connection type: 8-pin DIN Description: 2-cell conductivity/ATC electrode, 0.1 cell constant, SS sensor, epoxy shaft (700 and 2700 series) Diameter (metric) Body: 115 x 12 mm Length (metric) Cable: 1 m Material: Stainless steel rings with epoxy shaft Measuring ranges: lt;200 mu;S Type of filling solution: Ag/AgCl, refillable Unit Size: Each

Additional Information

No.	LD-2962
Manufacturer (Brand)	Thermo Scientific

