



Thermo Scientific Orion Lead Electrodes

Product Images

Description

Orion Lead Electrodes The Thermo Scientific Orion 9482BN Lead Half Cell and Orion 9682BNWP SureFlow Lead Electrode from the Ionplus Series can be used to measure lead ions in aqueous solutions quickly, easily, accurately and inexpensively. Connection types: Waterproof BNC Description: With waterproof BNC connector Deviation: Drift-free potentials Diaphragm: Sure-Flow™ free-flow diaphragm provides stability For use with (devices): ISE meters Type: Ion-selective lead electrode, waterproof BNC connector Unit Size: Each

Additional Information

No.	LD-3038
Manufacturer (Brand)	Thermo Scientific

