



Thermo Scientific Orion Lead Electrodes

Product Images

Description

Orion Lead Electrodes The Thermo Scientific Orion 9482BN Lead Half Cell and Orion 9682BNWP SureFlow Lead Electrode from the Ionplus Series can be used to measure lead ions in aqueous solutions quickly, easily, accurately and inexpensively. Connection types: BNC connector Description: With BNC connector Deviation: Drift-free potentials Diaphragm: Sure-Flow™ free-flow diaphragm provides stability For use with (devices): ISE meters Type: Ion-selective lead electrode, BNC connector Unit Size: Each

Additional Information

No.	LD-2996
Manufacturer (Brand)	Thermo Scientific

