



Thermo Scientific Orion Cadmium Electrodes

Product Images

Description

Orion Cadmium Electrodes The Thermo Scientific Orion 9448BN Cadmium Half Cell and Orion 9648BNWP Sure-Flow Cadmium Combination Electrode provide fast, easy, accurate and cost-effective measurement of cadmium ions in solution. Connection Types: Waterproof BNC Description: Model 9448BN ; Orion Cadmium Combination Diaphragm: Sure-Flow free-flow diaphragm provides stable, drift-free potentials For use with (devices): ISE meters Type: Ion selective combination electrode Unit Size: Each

Additional Information

No.	LD-3032
Manufacturer (Brand)	Thermo Scientific

