



Thermo Scientific Orion Calcium Electrodes

Product Images



Description

Orion Calcium Electrodes The Thermo Scientific Orion 9320BN and Orion 9720BNWP Ionplus Calcium Electrodes allow calcium ions in solution to be measured easily, accurately and inexpensively. Connection types: BNC connector Description: Combination Calcium Electrode with Waterproof BNC Connector Diaphragm: Free-flowing Sure-Flow junction assures stable, drift-free potentials For use with (instruments): ISE meters Type: Ion selective combination electrode Unit Size: Each

Additional Information

No.	LD-3047
Manufacturer (Brand)	Thermo Scientific

