



Thermo Scientific Orion Calcium Electrodes

Product Images



Description

Orion Calcium Electrodes The Thermo Scientific Orion 9320BN and Orion 9720BNWP Ionplus Calcium Electrodes make measuring calcium ions in solution easy, accurate and cost-effective. Description: Replacement Module for Half-Cell Calcium Electrode 9320BN Diaphragm: Free-flowing Sure-Flow junction assures stable, drift-free potentials For use with (instruments): ISE meters Type: Ion selective combination electrode Unit Size: Each

Additional Information

No.	LD-2958
Manufacturer (Brand)	Thermo Scientific

