



Thermo Scientific Orion Cyanide Electrodes

Product Images

Description

Orion Cyanide Electrodes The Thermo Scientific Orion 9406BN Cyanide Half Cell and Orion 9606BNWP Sure-Flow Cyanide Combination Electrode provide fast, easy, accurate and cost-effective measurement of cyanide ions in aqueous solutions. Connection types: Waterproof BNC Description: Model 9406BN ; Cyanide combination Diaphragm: Free-flowing liquid junction assures stable, drift-free potentials For use with (devices): ISE meters Type: Ion selective electrode Unit Size: Each

Additional Information

No.	LD-3024
Manufacturer (Brand)	Thermo Scientific

