



Thermo Scientific Orion Potassium Electrodes

Product Images

Description

Orion Potassium Electrodes The Thermo Scientific Orion 9319BN and Orion 9719BNWP Ionplus Potassium Electrodes make the determination of potassium in aqueous solution simple, accurate and cost-effective. Description: Replacement module for 9319BN half-cell electrode for potassium determination Deviation: drift-free potentials Diaphragm: Sure-Flow free-flow diaphragm provides stability For use with (instruments): ISE meters Type: Ion selective combination electrode Unit Size: Each

Additional Information

No.	LD-2955
Manufacturer (Brand)	Thermo Scientific

