



Thermo Scientific Orion Potassium Electrodes

Product Images

Description

Orion Potassium Electrodes The Thermo Scientific Orion 9319BN and Orion 9719BNWP Ionplus potassium electrodes make it easy, accurate and cost-effective to determine potassium in aqueous solution. Connection types: BNC connector Description: Sure-Flow combination electrode for potassium with waterproof BNC connector Deviation: drift-free potentials Diaphragm: Sure-Flow free-flow diaphragm provides stability For use with (instruments): ISE meters Type: Ion selective combination electrode Unit Size: Each

Additional Information

No.	LD-3045
Manufacturer (Brand)	Thermo Scientific

