



# Thermo Scientific Orion Potassium Electrodes

## Product Images

## Description

Orion Potassium Electrodes The Thermo Scientific Orion 9319BN and Orion 9719BNWP Ionplus potassium electrodes make the determination of potassium in aqueous solution simple, accurate and cost-effective. Description: Replacement module for Sure-Flow combination electrode 9719BNWP for potassium determination Deviation: drift-free potentials Diaphragm: Sure-Flow free-flow diaphragm provides stability For use with (instruments): ISE meters Type: Ion selective combination electrode Unit Size: Each

## Additional Information

|                      |                   |
|----------------------|-------------------|
| No.                  | LD-3044           |
| Manufacturer (Brand) | Thermo Scientific |

