



Thermo Scientific Orion Copper Electrodes

Product Images

Description

Orion Copper Electrodes The Thermo Scientific Orion 9429BN Copper Half Cell and Orion 99629BNWP Copper Electrode from the Ionplus Series can be used to measure copper ions in solution quickly, easily, accurately and inexpensively. Connection types: Waterproof BNC Description: With waterproof BNC connector Diaphragm: Free-flow diaphragm provides stable, drift-free potentials For use with (devices): ISE meters Type: Ion selective combination electrode Unit Size: Each

Additional Information

No.	LD-3030
Manufacturer (Brand)	Thermo Scientific

