



Thermo Scientific Orion Sodium Electrodes

Product Images

Description

Orion Sodium Electrodes The Thermo Scientific Orion Sodium Electrodes are designed for the rapid, accurate, and cost-effective measurement of sodium in aqueous solution. Connection types: BNC, waterproof Description: ROSS sodium combination electrode Type of filling solution: 900010 or 900012 for low Na+ content, ROSS filling solution for sodium electrodes: 50 ml Length (metric) Cable: 1 m Material: sodium Operating range (metric): 1 x 10-6 M to saturation 0.02 ppm to saturation For use with (instruments): ISE meters Operating temperature (metric): 0 ° to 100 °C Type: Ion selective combination electrode Unit Size: Each

Additional Information

No.	LD-2659
Manufacturer (Brand)	Thermo Scientific

