



Thermo Scientific Orion Nitrate Electrodes

Product Images

Description

Orion Nitrate Electrodes The Thermo Scientific Orion 9307BNWP and Orion 9707BNWP Ionplus Nitrate Electrodes provide simple, accurate and cost-effective measurement of nitrate in aqueous solutions. Description: Replacement module for 9307BNWP half-cell electrode for nitrate determination, 1/pack Deviation: drift-free potentials Diaphragm: Sure-Flow™ free-flow diaphragm provides stability For use with (instruments): ISE meters Type: Ion selective combination electrode Unit Size: Each

Additional Information

No.	LD-2949
Manufacturer (Brand)	Thermo Scientific

