



Thermo Scientific Orion REDOX/ORP/Temperature Electrodes

Product Images

Description

Orion REDOX/ORP/Temperature Electrodes The Thermo Scientific Orion combined redox/ORP/temperature electrodes are designed to measure REDOX potential and temperature and include a platinum REDOX sensor, an internal Ag/AgCl reference system, and a built-in thermistor. Description: low maintenance epoxy triode Material: epoxy For use with (instruments): Orion A+ Series Instruments Type: ORP/Temperature Electrode Unit Size: Each

Additional Information

No.	LD-2912
Manufacturer (Brand)	Thermo Scientific

