



Thermo Scientific Orion™ Bromide Electrodes

Product Images

Description

Orion™ Bromide Electrodes With Thermo Scientific™ Orion™ Bromide Electrodes, you can measure bromide ions in aqueous solutions quickly, easily, accurately, and cost-effectively. Connection types: Waterproof BNC Description: Model 9635 ; Orion Bromide Combination Diaphragm: Free-flowing Sure-Flow junction assures stable, drift-free potentials For use with (instruments): ISE meters Type: Ion selective combination electrode Unit Size: Each

Additional Information

No.	LD-3031
Manufacturer (Brand)	Thermo Scientific

