



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Amarell® immersion sensor Pt100 class A,  
stainless steel, 105x3,0mm, PVC cable approx.  
1m, for**

**€45.61** As low as  
**€39.49**  
**plus VAT &  
Shipping**

## Product Images

---



## Description

---

Immersion sensor Pt100 class A, stainless steel, 105x3,0mm, PVC cable approx. 1m, for ad 20 th

## Additional Information

---

No.	LL-0648
Manufacturer (Brand)	Amarell
VGKL number	681042010

