



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Analytik Jena® Sample holder P20 (aluminum, black)

**€212.00**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

Specimen holder P20 (aluminum, black) : For use with the SpeedMill PLUS\* and holding up to 20 Lysis Tubes, autoclavable, storage at up to -80 °C\*CAUTION! For SpeedMill PLUS up to S/N 10-3280 1250 a new closure is necessary for the use of the specimen holder P20, please order alternatively the "Set specimen holder P20 (aluminum, black)" (845-60055-0).Scope of deliverySpecimen holder P20 (aluminum, black)

## Additional Information

---

No.	AJ-0152
Manufacturer (Brand)	Analytik Jena
Transport temperature	Room temperature

