



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Analytik Jena® Set Sample holder P20 (aluminum, black)

€602.00
plus VAT &
Shipping

Product Images

Description

Set Specimen holder P20 (aluminum, black) : For use with the SpeedMill PLUS and holding up to 20 Lysis Tubes, autoclavable, storage at up to -80 °C, including new impact plates and closure* for the SpeedMill PLUS* ATTENTION: From S/N 10-3280 1250 the SpeedMill PLUS is equipped with new impact plates and closure as standard, alternatively the specimen holder can be ordered separately (845-60053-0).Scope of deliverySpecimen holder P20 (aluminum, black)Closure plateBase plate

Additional Information

No.	AJ-0149
Manufacturer (Brand)	Analytik Jena
Transport temperature	Room temperature

