



## Bandelin® SONOREX D 170 Lid

### Product Images

---



## Description

A suitable stainless steel lid for an ultrasonic bath protects the bath liquid from external contamination and prevents odor nuisance caused by the added cleaning or disinfection preparations. In addition, the noise emission during operation of the ultrasonic cleaner, which may be perceived as annoying, is reduced. The lids are designed in such a way that any condensation water produced is directed back into the tank.

## Additional Information

No.	BL-0063
Manufacturer (Brand)	Bandelin
VGKL number	960733006

