



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bandelin® SONOREX TECHNIK WD 300 Lid

**€483.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

A suitable lid for an ultrasonic bath protects the cleaning liquid from contamination and prevents odor nuisance from the cleaning or disinfecting liquid. In addition, the noise emission during operation of the ultrasonic cleaner, which may be perceived as annoying, is reduced. The lids are designed in such a way that any condensation water produced is directed back into the tank. Suitable for W 300

## Additional Information

No.	BL-0689
Manufacturer (Brand)	Bandelin
Transport temperature	Room temperature

