




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bandelin® SONOREX ZF 10 laboratory flask holder 8 pieces

**€154.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

For fixing laboratory flasks in hanging baskets Set consisting of 5 pcs. laboratory flask holders à 155 mm and 3 pcs. laboratory flask holders à 215 mm; stainless steel In order to carry out smooth sample preparation, appropriate conditions must be created for laboratory work. Samples should be homogenized, extracted or degassed quickly and reliably for subsequent analysis in laboratory vessels of various shapes and sizes. For this challenge, BANDELIN has developed practical holders for fixing laboratory flasks in suspension baskets for SONOREX ultrasonic baths. The holders guarantee the stability of the laboratory flasks during sonication. At the same time, they ensure optimum coupling of the samples to the contact liquid and thus low-loss ultrasound transmission into the sample to be sonicated. The laboratory flask holders can be mounted quickly and easily as longitudinal and transverse dividers at any position of the basket rim by means of small hooks. The compartment size can be individually adjusted according to the flask shape and size (10 - 250 ml) and can also be changed again in the blink of an eye if required. Without any effort or additional tools, the laboratory flasks can be inserted with an exact fit, so that sample preparation with ultrasound can be carried out safely and effectively. To ensure that the SONOREX ultrasonic bath continues to be filled at the optimum level and to avoid flooding of the laboratory flasks, we recommend the GV 10 handle adjustment.

## Additional Information

---

No.	BL-0217
Manufacturer (Brand)	Bandelin
Transport temperature	Room temperature

