



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bandelin® SONOPULS TS 419 Titanium sonotrode

€500.00
plus VAT &
Shipping

Product Images



Description

Beam diameter 19 mm, suitable for SH 400 G (HD 4400). Sonotrodes transfer the mechanical vibrations into the specimen, are thermostable, autoclavable and resistant to corrosive media. Sample volume, vessel diameter and desired amplitude determine instrument selection and sonotrode type. A high amplitude means particularly intensive sonication. The sonotrode is connected to step or booster horns. All BANDELIN sonotrodes are made of the titanium alloy TiAl6V4 (3.7165).

Additional Information

No.	BL-0282
Manufacturer (Brand)	Bandelin
Transport temperature	Room temperature

