



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bandelin® SONOREX TECHNIK FA 320 Filtration

**€5,032.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The FA filter unit is used to filter out dirt particles from the cleaning solution and enables the service life of the ultrasonic bath to be extended from a volume of 13 liters. The cleaning system is suitable for filtering aqueous cleaning media up to a maximum liquid temperature of 70 °C. The FA filter unit consists of two filter housings made of polypropylene (PP) and a magnetically coupled centrifugal pump, which are mounted together on a console made of stainless steel sheet. The pre-filter serves to protect the pump from coarse contamination. The fine filter is used for depth filtration. Further equipment includes a pressure switch, a ball valve, a switch box with ON button, a main switch and 2 meters of connection cable with safety plug. - The pressure switch switches off the pump when the fine filter is clogged. The - fine filter can then be replaced. - The ball valve can be used to adjust the flow rate. - The pump of the filter unit must be below the liquid level of the tub. Mains connection: 230 V EU plug

## Additional Information

---

No.	BL-0793
Manufacturer (Brand)	Bandelin
Transport temperature	Room temperature

