



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bandelin® SONOREX TECHNIK WG 180-2 Tub frame

€7,034.00
plus VAT &
Shipping

Product Images



Description

With the WG tub rack for 2 tubs, the lifting device with MB oscillation, which is available separately, can be moved between the cleaning stages. The cleaning baskets can thus be conveniently transported between the cleaning stages. The tub rack is used to realize multi-stage cleaning lines consisting of cleaning and rinsing basins as well as dryers.

Additional Information

No.	BL-0768
Manufacturer (Brand)	Bandelin
Transport temperature	Room temperature

