



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bohlender® SICCO desiccant cylinder 108 x 150 mm

€211.35
plus VAT &
Shipping

Product Images



Description

SICCO desiccant cylinder Aluminum cylinder with encapsulated silica gel and saturation indicator by color change from blue to pink on the top. Dimensions: Outer dimensions $a\varnothing \times h$ in mm 108 x 136 Filling capacity approx. g 750 Weight approx. g 980 suitable for

Additional Information

| | |
|-----------------------|------------------|
| No. | ER-2347 |
| Manufacturer (Brand) | Bohlender |
| Transport temperature | Room temperature |

