



Bürkle® Novartos Multi, V4A,100cm,3-fold,21 chambers,o.Ein.

Product Images

Description

Novartos Uno sample collector with two or three collection chambers, for one zone sampling. Novartos the pharmaceutical zone sampler for sampling according to unit-dose sampling FDA guideline. Double or triple sampling is now possible without multiple piercing of the lance into the bulk material, thus segregation is largely avoided. A decisive advantage of Novartos is the ease of cleaning. The complete filling and emptying of the sample chambers has been decisively improved by computer-optimized geometry; there are no gaps or dead spaces in which sample residues can accumulate. Novartos Uno is used for targeted sampling The precisely fitting sample inserts are easy to replace Inserts are available in volumes from 0 to 2.0 ml (Novartos Multi) or 0 to 4.0 ml (Novartos Uno) (alternative sample volumes at intervals of 0.1 ml each are available on request) The insertion depth is easy to read on the laser-engraved depth scale Stainless steel V4A (1.4404 / 1.4571) Scratch-free and electropolished 25 mm Ø Novartos was developed in cooperation with a multinational pharmaceutical/chemical company in Basel/Switzerland.

Additional Information

No.	KL-2308
Manufacturer (Brand)	Bürkle

