



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Carbolite-Gero® Base Inert Gas Package FHA/C 13/50/500

€3,570.00
plus VAT &
Shipping

Product Images



Description

The basic inert gas package enables working under inert gas and rough vacuum. The basic inert gas package consists in detail of the following items: Two gas-tight stainless steel flanges with quick-release fasteners and NW25 adapters
Two fiber insulation plugs with tubes up to 1400°C
Two flange holding sets mounted on the furnace housing
One long ceramic tube 610 (Mullite / Pythagoras) $\varnothing 46 \times \varnothing 40 \times 1225 \text{ mm}$ - Mullite (Pythagoras) = 60% Al₂O₃, 40% SiO₂ - Gas-tight - Tmax: 1500°C - Maximum thermal shock resistance (heating rate) 300°C/h (5°C/min)

Additional Information

No.	CX-0267
Manufacturer (Brand)	Carbolite-Gero
Transport temperature	Room temperature

