



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Carbolite-Gero® Base Inert Gas Package HTRH/-3 __/100/600

€10,110.00
plus VAT &
Shipping

Product Images



Description

The basic inert gas package enables working under inert gas and rough vacuum. The basic inert gas package consists in detail of the following items: Two gas-tight stainless steel flanges with quick-release fastener and NW40 adapter
Two fiber insulation plugs with tubes up to 1800°C
Two insulation sleeves to protect the working tube
Two table-based tube supports for horizontal operation
One long ceramic tube 799 (RCA / Alsint) $\varnothing 100 \times \varnothing 88 \times 1415 \text{ mm}$ - RCA (Alsint) = 99.7% Al₂O₃ - Gas-tight - Tmax: 1800°C - Two insulating sleeves to protect the working tube - Tmax under vacuum: 1450°C - Maximum thermal shock resistance (heating rate) 300°C/h (5°C/min)

Additional Information

No.	CX-0282
Manufacturer (Brand)	Carbolite-Gero
Transport temperature	Room temperature

