



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Carbolite-Gero® Base Inert Gas Package HTRV \_\_/70/500

**€7,660.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The basic inert gas package enables working under inert gas and rough vacuum. The basic inert gas package consists in detail of the following items: Two gas-tight stainless steel flanges with quick-release fastener and NW40 adapter  
Two fiber insulation plugs with tubes up to 1800°C Two insulating sleeves to protect the working tube Two tube holders for vertical operation One long ceramic tube 799 (RCA / Alsint)  $\varnothing 70 \times \varnothing 60 \times 1290 \text{mm}$  - RCA (Alsint) = 99.7% Al<sub>2</sub>O<sub>3</sub> - Gas-tight - Tmax: 1800°C - Two insulating sleeves to protect the working tube - Tmax under vacuum: 1450°C - Maximum thermal shock resistance (heating rate) 300°C/h (5°C/min)

## Additional Information

---

No.	CX-0288
Manufacturer (Brand)	Carbolite-Gero
Transport temperature	Room temperature

