



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Carbolite-Gero® Base Inert Gas Package HTRV- A __/70/120

**Price on
request
plus VAT &
Shipping**

Product Images



Description

The basic inert gas package enables working under inert gas and rough vacuum. The basic inert gas package consists in detail of the following items: Two gas-tight stainless steel flanges with quick-release fastener and NW40 adapter
Two fiber insulation plugs with tubes up to 1800°C Two insulating sleeves to protect the working tube Two tube holders for vertical operation One long ceramic tube 799 (RCA / Alsint) $\varnothing 70 \times \varnothing 60 \times 910 \text{mm}$ - RCA (Alsint) = 99.7% Al₂O₃ - Gas-tight - Tmax: 1800°C - Two insulating sleeves to protect the working tube - Tmax under vacuum: 1450°C - Maximum thermal shock resistance (heating rate) 300°C/h (5°C/min)

Additional Information

| | |
|-----------------------|------------------|
| No. | CX-0291 |
| Manufacturer (Brand) | Carbolite-Gero |
| Transport temperature | Room temperature |

