



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**Carbolite-Gero® CWF 11/13 basic controller E301  
chamber furnace 1100deg maximum  
temperature 13 l**

**€5,420.00  
plus VAT &  
Shipping**



## Description

---

Chamber oven 1100°C maximum temperature 13 l chamber volume 200x200x325mm (HxWxD) chamber size 655x435x610mm (HxWxD) external dimensions connection: 230 V, 50 Hz, 3,1 kW The picture shows an article of the same type. Size and equipment may differ from this type. PID controller One ramp to setpoint Timer function with 99 hours, programmable as delayed start or countdown timer. Optional CC-T1 touch screen controller (1 zone) (00050670) Overtemperature protection digital (OPTOTCF)

## Additional Information

---

No.	CX-0019
Manufacturer (Brand)	Carbolite-Gero
VGKL number	427111013
Transport temperature	Room temperature

