



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Carbolite-Gero® CWF 11/5 basic controller E301
chamber furnace 1100degrees maximum
temperature 4.7 l**

**€4,510.00
plus VAT &
Shipping**



Description

Chamber oven 1100°C maximum temperature 4,7 l Chamber volume 135x140x250mm (HxWxD) Chamber size 585x375x485mm (HxWxD) Outer dimensions Connection: 230 V, 50 Hz, 2,4 kW The picture shows an article of the same type. Size and equipment may differ from this type. PID controller One ramp to setpoint Timer function with 99 hours, programmable as delayed start or countdown timer. Optional CC-T1 touch screen controller (1 zone) (00050670) Digital overtemperature protection (OPTOTCF)

Additional Information

No.	CX-0021
Manufacturer (Brand)	Carbolite-Gero
VGKL number	427111005
Transport temperature	Room temperature

