



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Carbolite-Gero® TSO1 11/600 1100degrees, 600mm EPC3016P1 PID controller TSO1 11/600 1100degrees, 600mm

**€15,140.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

TSO1 11/600 1100°C, 600mm Heated length, single zone rotary reactor Tube furnace (220 - 240 V single phase)  
EPC3016PI with 24 free programmable segments (e.g. 12 ramps and 12 hold times) - Self-optimization - Thermocouple  
break protection - Electronic setting limit - Digital temperature setting - Digital actual and set point display - Max. 2  
control tracks (e.g. for optional solenoid valves) - Ethernet connection on the rear panel of the control housing  
Standard TSO/600 quartz vessel with smooth inner surface Optional surcharge for CC-T1 of EPC3016PI (00053272)  
Digital overtemperature protection (00053325) Upgrade TSO/600 to quartz vessel with fluted inner surface

## Additional Information

---

|                       |                  |
|-----------------------|------------------|
| No.                   | CX-0406          |
| Manufacturer (Brand)  | Carbolite-Gero   |
| Transport temperature | Room temperature |

